

Protégé

(Protégé)

Introducing Aurelia® Protégé™ stretch nitrile exam gloves. **The comfort and fit of latex, without the possible allergic reaction.** It's so stretchy that your long hours seem to fly by. And you won't believe the sensitivity from such a durable, yet lightweight glove. With the new, ergonomic Protégé, you can work smarter not harder.



NEW!



**Improved Protection
Amazing Comfort**



AURELIA®
Sign of Comfort™



www.aureliagloves.com.sg | Tel: 65 6911 3208 | Fax: 65 6911 3203

Protégé

(Protégé)

Product Description

Aurelia® Protégé™ stretch nitrile powder free examination gloves are made from 100% nitrile, ambidextrous, sky blue color and non-sterile. Each batch of gloves is air-tested to detect pinholes. Aurelia gloves conform and are tested to the highest international standards. For more information, visit www.aureliagloves.com.sg

- 100% Latex Free, non-latex allergenic (Worry Free, Latex Free)
- Enhanced nitrile formulation rivals latex comfort and ergonomics
- Fingertip-textured finish for precise gripping
- Greater chemical resistance compared to other material
- Ambidextrous
- Non-Sterile, Single-Use
- 100 Gloves per Box*



Aurelia® Protégé™ Specifications

TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
Small	32.6	725
Medium	36.1	670
Large	34.9	660
Extra Large	35.8	680
ASTM D6319-00a	14 MPa minimum	500% minimum

AGED TENSILE TEST RESULTS (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)
Small	30.8	620
Medium	34.2	635
Large	33.6	620
Extra Large	33.4	610
ASTM D6319-00a	14 MPa minimum	400% minimum

DIMENSIONS

	Length (inches)	Width (inches)	Thickness (Single-wall) Palm (Mils)	Finger (Mils)
Small Aurelia Protégé ASTM D6319-00a	9.5 minimum 9.1 minimum	3.5 ± 0.08 3.2 ± 0.4	4.0 ± 1.2 2.0 minimum	5.2 ± 1.2 2.0 minimum
Medium Aurelia Protégé ASTM D6319-00a	9.5 minimum 9.1 minimum	3.7 ± 0.08 3.7 ± 0.4	4.0 ± 1.2 2.0 minimum	5.2 ± 1.2 2.0 minimum
Large Aurelia Protégé ASTM D6319-00a	9.5 minimum 9.1 minimum	4.1 ± 0.08 4.1 ± 0.4	4.0 ± 1.2 2.0 minimum	5.2 ± 1.2 2.0 minimum
Extra Large Aurelia Protégé ASTM D6319-00a	9.5 minimum 9.1 minimum	4.5 ± 0.08 4.5 ± 0.4	4.0 ± 1.2 2.0 minimum	5.2 ± 1.2 2.0 minimum

Ordering Information

To order from your preferred distributor, specify size and item number. Packing: 100 pieces per dispenser box*, 10 dispenser boxes per shipping carton.

*Note: Extra-Large 90 gloves by weight per carton

ITEM NUMBER

Small	Medium	Large	Extra Large
93996	93997	93998	93999

Authorized Aurelia Distributor:

