


Sonic®

Aurelia® Sonic Nitrile® Indigo Blue
Nitrile Examination Gloves



- New ultra comfort formulation
- 2.2 mil thickness for improved tactile sensitivity
- Finger textured finish
- Cost saving & Space saving: 100 gloves per box
- Soothing indigo blue color



AURELIA®
Sign of Comfort™

www.aureliagloves.com.sg | Tel: 65 6911 3208

Sonic®

Aurelia is proud to introduce our latest innovation in personal protection, Aurelia Sonic® Nitrile Gloves. Featuring our revolutionary ultra-thin formulation Sonic® will adapt so seamlessly to your hands that you'll be working faster than the speed of sound! Each box of Aurelia Sonic® contains 100 machine packed soothing indigo blue colored gloves. We guarantee Sonic® will be the next big BOOM for healthcare professionals! For more information visit www.aureliagloves.com.sg



Aurelia® Sonic® Specifications

TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	30.0	680
Small	31.1	680
Medium	31.0	680
Large	31.2	680
Extra Large	31.3	680
ASTM D6319	14 MPa minimum	500% minimum

AGED TENSILE TEST RESULTS (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	30.0	620
Small	30.0	620
Medium	30.0	620
Large	30.0	620
Extra Large	30.0	620
ASTM D6319	14 MPa minimum	620% minimum

DIMENSIONS

	Length (inches/mm)	Thickness (Single-wall) Palm (mils/mm)	Finger (mils/mm)
Extra Small Aurelia Sonic ASTM D6319	9.4 / 240 8.7 / 220 minimum	2.2 / 0.06 2.0 / 0.05 minimum	3.2 / 0.08 2.0 / 0.05 minimum
Small Aurelia Sonic ASTM D6319	9.4 / 240 8.7 / 220 minimum	2.2 / 0.06 2.0 / 0.05 minimum	3.2 / 0.08 2.0 / 0.05 minimum
Medium Aurelia Sonic ASTM D6319	9.4 / 240 8.7 / 220 minimum	2.2 / 0.06 2.0 / 0.05 minimum	3.2 / 0.08 2.0 / 0.05 minimum
Large Aurelia Sonic ASTM D6319	9.4 / 240 8.7 / 220 minimum	2.2 / 0.06 2.0 / 0.05 minimum	3.2 / 0.08 2.0 / 0.05 minimum
Extra Large Aurelia Sonic ASTM D6319	9.4 / 240 8.7 / 220 minimum	2.2 / 0.06 2.0 / 0.05 minimum	3.2 / 0.08 2.0 / 0.05 minimum

Ordering Information

To order from your preferred distributor, specify size and item number. Packing: 100 pieces per dispenser box, 10 dispenser boxes per shipping carton.

ITEM NUMBER

Extra Small	Small	Medium	Large	Extra Large
9377B5	9377B6	9377B7	9377B8	9377B9

Authorized Aurelia Distributor:



AURELIA
Sign of Comfort™

MAXTER Healthcare PTE LTD

One Commonwealth Lane, Unit# 07-01 and 07-02, Singapore 149544 | Tel: 65 6911 3208 | Fax: 65 6911 3203
www.aureliagloves.com.sg | Email: info@aureliagloves.com.sg