


Transform[®]

Treat your hands like superstars with Aurelia's newest comfort and performance series nitrile gloves. Introducing Transform[®] - our newest nitrile gloves that combine all the best benefits of a nitrile glove in a proprietary formulation that take strength, comfort and performance to an all new level.



Perform Like A Superstar.
Transform Your Game.

- Ultra-thin for superior tactile sensitivity
- Combines long-lasting comfort, durability and decreased hand fatigue
- Improved barrier protection
- New TransBlue[®] color



AURELIA[®]
Sign of Comfort[™]

www.aureliagloves.com.sg | Tel: 65 6911 3208

Transform[®]

Product Description

Treat your hands like superstars with Aurelia's newest comfort and performance series nitrile gloves. Introducing Transform[®] - our newest nitrile gloves that combine all the best benefits of a nitrile glove in a proprietary formulation that take strength, comfort and performance to an all new level.



Aurelia[®] Transform[®] Specifications

TENSILE STRENGTH RESULTS (before aging)

Tensile Strength (MPa)	Elongation at break (%)
min 14	min 500

TENSILE STRENGTH RESULTS (after aging)

Tensile Strength (MPa)	Elongation at break (%)
min 14	min 400

DIMENSIONS

Length (mm/in)	Cuff Thickness (mils)	Palm Thickness (mils)	Finger Thickness (mils)
240/9.5	2	2.8	4

Aurelia[®] Transform[®]

Material: Nitrile (Latex Free)

Glove Design: Ambidextrous

Glove Feature: Powder Free and Finger Textured

Color: Trans Blue

Packaging: 2000 gloves/case

10 boxes 200 gloves/case



Ordering Information

To order from your preferred distributor, specify size and item number. Packing: 200 pieces per dispenser box, 10 dispenser boxes per shipping carton.

ITEM NUMBER

Extra Small	Small	Medium	Large	Extra Large
98895	98896	98897	98898	98899

Authorized Aurelia Distributor:



AURELIA[®]
Sign of Comfort[™]

MAXTER Healthcare PTE LTD

One Commonwealth Lane, Unit# 07-01 and 07-02, Singapore 149544 | Tel: 65 6911 3208 | Fax: 65 6911 3203
www.aureliagloves.com.sg | Email: info@aureliagloves.com.sg