


Velocity

(Velocity)

Introducing Aurelia® Velocity™ triple washed, polymer-coated, powder free latex examination gloves. **Perhaps the finest glove available.** Advance technology for maximum purity. Maximum quality. And wait until you feel the grip that Velocity's polymer-coating delivers — wet or dry.



NEW AND IMPROVED!

**Maximum Purity
Superior Grip**



AURELIA®
Sign of Comfort™



www.aureliagloves.com.sg | Tel: 65 6911 3208 | Fax: 65 6911 3203


Velocity[®]

(Velocity)

Product Description

Aurelia[®] Velocity[™] polymer-coated, powder free latex examination gloves are made from 100% natural rubber, ambidextrous, creamy white color and non-sterile. Aurelia gloves conform and are tested to the highest international standards. For more information, visit www.aureliagloves.com.sg

- 100% visually inspected for maximum quality assurance and performance.
- Advanced triple washed process that significantly reduces risk of latex allergy, super low protein content.
- New dispensing-friendly package, with layered packing.
- Polymer-coating technology delivers excellent grip, enhances donning
- Powder Free
- Ambidextrous, Non-Sterile, Single-Use
- 100 Gloves per Box



Aurelia[®] Velocity[™] Specifications

TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	30.0	970
Small	31.1	970
Medium	31.0	970
Large	31.2	960
Extra Large	31.3	970
ASTM 3578:05	18 MPa minimum	650% minimum

AGED TENSILE TEST RESULTS (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	28.0	900
Small	28.0	920
Medium	27.5	910
Large	27.3	890
Extra Large	27.9	900
ASTM 3578:05	14 MPa minimum	500% minimum

DIMENSIONS

	Length (inches/mm)	Thickness (Single-wall)	
		Palm (mils/mm)	Finger (mils/mm)
Extra Small Aurelia Velocity ASTM 3578:05	9.4 / 240 8.7 / 220 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
Small Aurelia Velocity ASTM 3578:05	9.4 / 240 8.7 / 220 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
Medium Aurelia Velocity ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
Large Aurelia Velocity ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
Extra Large Aurelia Velocity ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum

Ordering Information

To order from your preferred distributor, specify size and item number. Packing: 100 pieces per dispenser box, 10 dispenser boxes per shipping carton.

ITEM NUMBER

Extra Small	Small	Medium	Large	Extra Large
28225	28226	28227	28228	28229



Authorized Aurelia Distributor: